Signia Motion Charge&Go X

Standard, Power and Super Power Benefits Overview



Provider Benefits

- Complete lineup of Rechargeable, Signia Xperience BTEs with the world's first rechargeable Super Power for added patient convenience.
- Smaller, sleek designs from previous models for added comfort.
- Advanced benefits with Power and Super
 Power options.
 - Upgrade to the Dry&Clean charger to keep hearing aids in top working conditions and reduce appointments.
 - **Option T-Coil** for patients leveraging looping systems for enhanced listening.

Wearer Benefits

- Comfort and convenience without limits technology to support moderate, severe and profound hearing losses.
- Enhanced situational awareness with adaptive technology that automatically adjusts to movements and changes around you.

• Compact chargers with long-lasting wear time to get up to 60+ hours of wear time on a single change. Charger models and options vary by model type.

Key Features*

- Up to 61 hours of use on a single charge with SP, 30 hours for P and 24 hours for Standard (based on a full 4-hour charge)
- 30 minute quick charge for an average of 5 hours of use
- Automatic on/off functionality with charger
- Extended frequency range offers wearers a competitive edge over other BTEs
- Ultra HD e2e
- Remote Care Ready

Compatible with:



StreamLine TV™



StreamLine Mic[™] miniPocket[™]



24/7 Digital

Assistant

7/5/3/2/1 and sDemo



Bluetooth

Connectivity

Ams Acoustic Motion





Processing



Sensors

Li-ion

Rechargeability

Available Colors:



Signia App is equipped with with Signia Face Mask Mode Improves audibility for those communicating with someone wearing a face mask - all with the touch of a button in the Signia App.

*Refer to Feature Overview Chart in the Signia Platform Brochure for more details

Туре	Motion Charge&Go	Motion Charge&Go Power	Motion Charge&Go Super Power
Available Performance Level	7/5/3/2/1	7/5/3/2/1	7/5/3/2/1
Matrix	60 dB/133 dB SPL 60 dB/125 dB SPL	77 dB/135 dB SPL 66 dB/130 dB SPL	82 dB/140 dB SPL
Battery Size	Lithium-ion power cell		
Features			
Frequency channels/ compression controls	7X - 48/20 5X - 32/16 3X - 24/12 2X - 16/8 1X - 16/8	7X - 48/20 5X - 32/16 3X - 24/12 2X - 16/8 1X - 16/8	7X - 34/16 5X - 32/16 3X - 24/12 2X - 16/8 1X - 16/8
Hearing Programs	7/5/3-6 2/1-4	7/5/3-6 2/1-4	7/5/3-6 2/1-4
Bluetooth [®] Streaming	•	•	•
TwinPhone ¹	٠	٥	•
TeleCare	٠	•	•
Tinnitus therapies ²	٠	٠	•
CROS/BiCROS ³	٠	٠	٠
Rocker Switch ⁴	٠	٠	٠
T-coil		٠	٠

Туре	Motion Charge&Go	Motion Charge&Go Power	Motion Charge&Go Super Power	
App/Accessory Compatibility				
Signia App	٠	•	٠	
StreamLine TV/ StreamLine Mic	٥	٠	٠	
miniPocket	•	•	•	
Noahlink Wireless	٠	•	٠	
Fitting set	Ear Wear 3.0	Ear Wear 3.0	Ear Wear 3.0	

All Motion Charge&Go X models work exclusively with the new Ear Wear 3.0. Contact your Signia Representative to learn more about the new Ear Wear 3.0. **Upgrade to the Dry&Clean Charger at purchase for Power and Super Power.**



Signia Motion Charge&Go X



Signia Motion

Charge&Go X P





Signia Motion Signi Charge&Go X SP Char f

Signia Dry&Clean Charger option for P & SP

(1) Available for levels 7 / 5 / 3 and activated via rocker switch, app or accessory

(2) Available for 7/5/3/2

(3) Compatible with CROS Pure Charge&Go X for levels 7/5/3

(4) Programmable via Connexx

The Bluetooth[®] word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. SIV-9847 SH202109 Copyright © 2021 Signia GmbH. All rights reserved.

